

No Paste. No Splatter. No Grit.

Experience all the benefits of the New Butler Paste-Free prophylaxis angle.

A few simple tips for optimal use.

1. The cup on the Paste Free prophylaxis angle may feel slightly stiff compared to other prophylaxis cups. This is by design and you will find that it will soften after just a few seconds of use. Start polishing in an upper quadrant where moisture from the saliva duct will begin to soften the cup, which will enhance the adaptation and the cleaning process.
2. The Paste-Free prophylaxis cup uses a medium grit pumice which means it is perfect for most patients. It may not be appropriate for patients with heavy staining.
3. Use normal pressure when cleaning with the new Paste-Free prophylaxis angle. Studies show it cleans as well as paste without all the splatter and grit. If you experience stalling or spinning, you may be pushing too hard on the teeth. Simply attach another angle and try again with normal pressure.



Item# 1210P